



LUMICEUTICALS

Shine With Harmonic Light



Harmonic Light

Light Pulsed at Harmonic Resonant Frequencies

Supports cellular repair & regeneration as
the ultimate biological NUTRIENT
&
TUNES the mind, body and spirit.



Non-invasive

Supports Self Healing



FDA Cleared For Pain Relief, Circulation
Enhancement, Fine Lines & Wrinkles and Injury
Recovery

Encompasses All Energetic Wellness Modalities

LumiCeuticals – Ultimate Self Healing


Harmonic Light is an Optimal Environment that Supports the Body's Natural Ability to Heal, Regenerate & Thrive




*"What you seek
is seeking you"
Rumi*

Light Therapy:

- Used all over the world
- **Thousands** of hospitals, clinics and private practices
- **30+** years global use
- **2,500+** research studies of light technology on the human body
- Every medical discipline is using or researching light energy
- Research for yourself the latest findings at **pubmed.com**
- Also called *Low Level Light Therapy, Photobiomodulation*





Relieve Pain, Improve Circulation & Enhance Recovery FDA Over The Counter Cleared



Increased Circulation and the General Temporary Relief of Pain & the Associated Effects of:

- Muscle/Joint Aches & Stiffness
- Arthritis – Rheumatoid and Osteoarthritis
- Tendonitis Pain
- Muscle Spasms & Strains
- Back Pain
- Repetitive Use Injuries like Carpal Tunnel
- Increased Rate of Injury Recovery & Wound Healing
- Sports Injuries
- Neuropathy

- Fine Line & Wrinkle Treatment with our Facial Rejuvenator



The Body Is A Living Photocell

We Are Designed To Receive Light

1. Eyes - Reticular Activating System



2. Through the Skin - Transdermal
From deep brain and bone.



3. In Our Field - Quantum Effect.




Receiving through the Eyes

A single photon hitting the back of the retina lights up the brain.

We receive light through our eyes into our light meter, the pineal gland.

Light synchronizes us with our natural world and consciousness.




The Pineal Gland

- **Sleep Cycles**
- **Mood**
- **Hormones**
- **Intuition**
- **Memory**
- **Creativity**
- **Consciousness**



The Endocrine System



Nervous System Balance

**“95% of all dis-ease
is caused by stress.”**

-National Institute of Health, 2010



**Stress of any type takes our
body out of it's healing
response.**

**The ANS is controlled by
light and sound impulses.**

Harmonic Light is Meditative

“When pathways are clear, flowing and carrying light information factors, the internal parasympathetic system is in control.

Thus, light therapies are a priori forms of meditation & healing.”


Gary Buchanan -Sonatherapy



Reaching A Meditative State

* Effects of Meditation:

- On:** genes associated with energy metabolism
 mitochondrial function
 insulin secretion
 telomere maintenance
- Off:** those involved in inflammation





*These effects were more pronounced and consistent for long-term practitioners.

* **IMPLICATIONS:** People who practice simple meditation aren't "just relaxing," they're experiencing "a **specific genomic response that counteracts the harmful genomic effects of stress.**"

* **CONCLUSION:** "Do it for years," said researchers at the Benson-Henry Institute for Mind/Body Medicine at Mass Gen, "and then these effects are quite powerful in how they change your gene activity."

Transdermal Effect of Harmonic Light Impact Skin to Bone





Red: Studies Show it Supports Soft Tissue Health

Muscle (Heart), Fatty Tissue, Gums, Wounds, Elevated Mood, Skin, Organ, Brain Circulation, Nitric Oxide Induction, Inflammation, Pain

Near Infrared: Studies Show it Supports Soft and Hard Tissue Health


Bones, Joints, Tendons, Teeth, Nitric Oxide Induction: Pain, Inflammation, Lymphatics, Brain. Deeper Penetration Depth in the Body.

Blue: Studies Show It is Cleansing & Supports Skin & Liver Health


Collagen, Mood Calming, Inflammation, Nervous System Calming, Auto-Immune Support, Liver, Anti-Pathogenic Properties: Anti-Viral, Anti-Fungal, Anti-Bacterial

“The Function of our Entire Metabolism is Dependent on Light “

Albert Fritz Popp




Light & Energy Production



LIGHT = The Ultimate Biological Nutrient

- Light donates photons to cells that **need** energy
- Photons increase the production of ATP (150%)
- Adenosine Triphosphate, or ATP, is the energy of the cell
- ATP energy governs 80-90% of the cells' activities



Athletic Performance




Muscle Recovers Better With Light

Lasers Med Sci. 2016 Apr;31(3):531-8. doi: 10.1007/s10103-016-1888-9. Epub 2016 Feb 12. Light-emitting diode therapy increases collagen deposition during the repair process of skeletal muscle.

Light-emitting diode therapy in exercise-trained mice increases muscle performance, cytochrome c oxidase activity, ATP and cell proliferation. © 2014 WILEY-VCH Verlag GmbH & Co. KGaA, Weinheim. *J Biophotonics.* 2015 Sep;8(9):740-54. doi: 10.1002/jbio.201400087. Epub 2014 Nov 6.

Light & Oxygen




Oxygen



- * A cleaner that assists immune system, cell signaling and homeostasis – balance in the body.
- * Too much ROS, Reactive Oxygen, can be damaging
- * UV and excess Heat can create too much ROS.
- * Applying the right light to the body “triggers the production of reactive oxygen species... without inducing DNA damage”

* Article (PDF Available) in Journal of Biomedical Optics 19(4):48002 · April 2014

Light & Alkalinity



“Cancer grows in oxygen deprived acidic tissue.”


“Disease cannot thrive in an alkaline body.”

Dr Otto Heinrich Warburg

Winner of the 1931 Nobel prize in Physiology



Light & Circulation



What *is* Nitric Oxide?

A signaling molecule produced in the endothelial cell linings along 100,000 miles of blood vessels and capillaries in the human body.

The Nobel Prize was awarded in 1998 for understanding its roles.



Role Of Nitric Oxide

Decreases pain

Keeps blood vessels open, elastic and functioning properly

Lowers Blood Pressure

Increases Oxygen Levels Throughout The Body

Ensures Blood Flow To The Brain

Increases energy production

Stimulates collagen production


Defends Against Tumor Cells

Decreases Numbness

Increases Lymphatic Activity

Kills bacteria

Increases blood flow and prevents the build up of plaque in the arteries



Studies Have Found That Nitric Oxide

- * Protects against GI ulcers
- * Protects against Alzheimer's
- * Prevents male pattern baldness & E.D.
- * Prevents early stages of Atherosclerosis
- * Prevents vascular complications of Diabetes
- * Improves heart function whether slower or faster by carrying oxygen in and carbon dioxide out.
- * Stimulates body to release Human Growth Hormone (HGH)
- * Improves Immune System Function
- * Improves Cardiovascular Activity

"I have come to see Nitric Oxide as the physical expression of Prana – your life force. It is the "Wow" molecule that continually resets our ability to connect body, mind, heart, and spirit. To activate this, we need to care for ourselves with plenty of sleep, exercise, life-giving foods, letting go of resentment, and opening up to affirmation and love."

{Christiane Northrup, M.D.}

Scientific Studies on Nitric Oxide

- **Aspirin for pain!**

Aspirin induces nitric oxide

British Journal of Pharmacology 2004; 143: 159-165

- **Fibromyalgia Pain**

Fibromyalgia caused by Abnormal microcirculation

Rheumatology (Oxford) 2000; 39: 917-921

- **Nitric oxide regulates wound healing**

J Surg Res 1996; 63: 237-240

- **Reversal of diabetic peripheral neuropathy with phototherapy** Age Ageing, 2006; 35: 11-16

- **ADHD**

In addition, evidence strongly suggests a reduction in blood flow in the prefrontal cortex for ADHD individuals.

- **Nitric oxide inhibits HIV-1 replication in human astrocytoma cells**

Biochem Biophys Res Commun. 1999; 254: 200-2.

- **Modulation of tendon healing by nitric oxide**

Murrell GA, Szabo C, Hannafin JA, Jang D, Dolan MM, Deng XH, Murrell DF, Warren RF. Inflamm Res 1997; 46 19-27.

Supporting Wound Healing

Diabetic wounds.

Tried to heal them for 9 months.





Completely healed after

9 Weeks of Daily Use With Our System



Light & Regeneration





Light
=
Photons
=
Energy
=
Regeneration



REGENERATION IN THE BRAIN

Harvard Medical School – Cellular Health & Light

- **Traumatic Brain Injury: Low-level light ... for effective therapy of injured brain.**

Vascular damage occurs frequently at the injured brain. “*Strikingly, these adverse events were reversed significantly by noninvasive exposure of injured brain to low-level light (LLL).*”

2015 May 13. doi: 10.1038/jcbfm.2015.87.

- **Harvard Psychiatrists: The Psychiatry of Light.**

“Efficacy extends beyond seasonal affective disorder and includes non seasonal depression and sleep disorders, with emerging evidence for a role in treating attention-deficit/hyperactivity disorder, delirium, and dementia.”

[Harv Rev Psychiatry](#). 2015 May-Jun;23(3):188-94. doi: 10.1097/HRP.000000000000078.


Skin Rejuvenation



Studies Have Shown That Light

- * **Stimulates Collagen Production** = younger looking, more plump skin.
- * **Repairs, Renews, Regenerates Skin Cells**
- * **Reduces the Appearance of Wrinkles, Diminishes Scarring, Shrinks Pores, Reduces Discoloration, and Speeds up the Healing Time of Blemishes.**
- * Combining blue light with red light helps **kill bacteria** associated with acne.
- * ***To help stop or reverse sun damage, Dr. Sarnoff and The Skin Cancer Foundation recommend the following:***
 - “Try a light-emitting diode (LED) – A lower energy LED can help promote collagen production and reduce fine lines and wrinkles.”***

Light & Detox



Light energizes cell membranes so that nutrients can better nourish cells and help them thrive and toxins can be eliminated.

Weight Loss & Metabolism Research

- ***Weight Loss & Metabolism:***

Phototherapy to reduce body fat, insulin resistance and "metabolic inflexibility" related to obesity in women undergoing weight loss treatment. *Lasers Surg Med.* 2015 Oct;47(8):634-42. doi: 10.1002/lsm.22395. Epub 2015 Jul 29.

- ***Bright light for weight loss: results of a controlled crossover trial.***

Morning bright light treatment reduces body fat and appetite in overweight women and may be included in weight control programs.

PMID: 23429094 [PubMed - indexed for MEDLINE] *Obes Facts.* 2013;6(1):28-38. doi: 10.1159/000348549. Epub








Pulsing Light at Harmonic Resonant Frequencies

Pulsed light is more penetrating and efficacious. Absorption and activation of the light take place much faster.” Gary Buchanan Sonotherapy



The body is held together by sound. The presence of disease indicates that some sounds have gone out of tune.

Deepak Chopra



MOZART SYMPHONY

IMAGINE - JOHN LENNON

LOVE

FUJIWARA BEFORE PRAYER



PEACE

THANK YOU

I WILL KILL YOU

FUJIWARA AFTER PRAYER

Pulsing Light at Harmonic Frequencies

“Sickness results when cells, molecules or particles of matter are out of their normal resonance or vibratory pattern. By repeatedly exposing damaged tissue to the normal resonance frequencies associated with that tissue, healing often occurs, sometimes quite rapidly.”



Dr. Paul Nogier
French Neurologist



Nogier Discovered Frequencies For Tuning The Physical Body


1. 73hz	2. 147hz	3. 294hz	4. 587 hz	5. 1174 hz	6. 2349 hz	7. 4698 hz
Cells, hormones circulation Immunity Mood Bones	inflammation, infection, scar tissue Mental health	nerves, wounds, skin, eyes, acute	GI tract, neuropathy, pancreas Metabolism Lymph Liver Circulation	bones, joints, ligament, tendon, muscles lungs	Brain, Relaxation Plateaus, Chronic	pain, nerves, spurs calcification

Every LUMI System Can Pulse Light at These Frequencies Which Correspond to Settings 1-7.

<p>A Activate Balanced Energy</p>	<p>B Be Calm De-Stress Emotional Balance Spiritual Tuning</p>	<p>C Rejuvenate Physical Tuning Elemental Tuning to Nature's Harmonics</p>
<p>Nogier Frequencies 1-7</p>	<p>Harmonies of 3,6,9 Solfeggio and Sine Frequencies</p>	<p>Harmonies of 432hz Pythagorean Golden Ratio 3:2</p>

Sweeping Settings A, B, C Deliver Pulsed Light in a Series of Frequencies During a 21 minute Relaxing Session





We are just beginning to prove the ancient wisdom that all life is organized according to divine geometry.

Reported Results

Reported Observations from The Therapeutic Learning Center Of Utah
586 Patients of Nancy Foster PhD, APRN Using our Light Therapy Systems Spring 2014

- * Typical patients experienced one 20 minute B setting sessions of light once a week during the study time period. A minimum of 4 sessions was typical.
- * 56 patients reported acceleration of the healing processes of wounds, burns, and skin issues like pustules.
- * 148 patients reported 50 to 75% reduction in pain.
- * 220 patients were able to discontinue their anti-depressant medication while under direct supervision.
- * 20 patients that experienced social anxiety reported less isolation to home and increased motivation and ability to engage with family and church.

- 10 patients experiencing **extreme anxiety and panic attacks**, refractory to medications **were prevented from having hospitalization.**
- 6 patients with ADHD and experiencing **learning problems reported increased focus, concentration, and memory. They also reported improved school grades and relationships.** No medications used
- Results: In 4 B setting sessions: 18 patients, whose lab studies showed decreased Glutathione levels and associated debilitating fatigue **showed 50 % increase in Glutathione and 15 to 25% increase in reported energy levels.**
- Overall Dr. Foster and the licensed massage therapist that she has administer the light therapy, Ally Sanchez, found that our light systems gave them : **“the ability to work with individual patients, multiple health issues, with unprecedented ease, speed and precision.”**

Tuning to the Pythagorean 3:2 Ratio

1	3	9	27	81	243	729
2	6	18	54	162	486	
4	12	36	108	324	972	
8	24	72	216	648		
16	48	144	432	-flowing water, mala beads, bowls.		
32	96	288	864			
64	192	576				
128	384	1152				
256	768					

Our LUMI C Setting

Macro Minerals

Boron 16 hz
Magnesium 36hz
Sulfur 48hz
Potassium 54hz
Calcium 64hz
Vanadium 72 hz
Iron 81 Hz
Silver 162 Hz

Micro

Minerals:
Iridium 288hz
Germanium 108hz
Molybdenum 144
hz

Source: Dr. A.
True Ott, PhD



432hz

Programmable Frequencies

- * Lumi Vibe Pro6 controller can pulse light at 1Hz- 9999Hz.
- * Easily custom program your system to take advantage of Frequency Suggestions from The Electroherbalism Frequency Lists and noted Health Professionals.
- * Research is showing that harmonic pulsed light can reprogram cell gene expression, "activate natural repairing abilities, and counteract cellular aging processes, paving the way toward unprecedented strategies of regenerative medicine." Glob Adv Health Med. 2014 Mar; 3(2): 40–55.
- * 92 Frequency Self-Sweeping Programs at the Touch of A Button




LIGHT ENERGY CUSTOM SWEEP POSSIBILITIES


Compiled From the Electroherbalism Frequency Handbook and True Ott PhD

<u>Whole Body Support</u>	
11 General Power Sweep	432 787 880 1550 2720 5000 9999
12 Grounding Balance - Ott	180 360 432 540 720 900 1080
13 Abundance	343 2794 4644 4790 4903 4916 4952
14 Transformation	2664 3330 3996 2475 855 1584 1746
15 Chakra Balancing	228 606 728 512 768 288 432
16 Instant Vacation	120 432 5000 9999
17 Master Healing Support	4170
___ Chakra 2	912 606 183 128 192 144 432
___ Chakra Wand Sweep	194 210 126 136 141 221 172
___ DNA Repair Support	432 731 732 432
___ Energize	1056 2008 9999 424 27
___ Schumann Wave	10 14 20 26 33 39 45
___ Ultimate Health 432	432 27 54 108 216 432 1080 27 144 432
<u>Mood and Sleep Support</u>	
18 Peace & Calming	1009 4601
19 Insomnia Support	48 802 880 1550 6000 8778 9999
20 Depression /Anxiety Support	9999 5000 3176 800 787 73
___ Emotion Sweep	4992 7240 3727 4582 4743 4875 4929
___ Joy & Happiness	432 4604 4655 4634 4589
___ Sleep Apnea	15 10 20 35 40 80 125 150 160 240 320 635 763 2720
<u>Immune & Cleanse Support</u>	
21 A+ Immune Support	337
22 Allergies	146 290 337 555 690 727 787
23 AutoImmune Support	328 880 728 5000 2720 2180 120
24 Colds/Flu Support	800 880 8700 7760 4888 5000 9999
25 Chem Detox	664 7344 2842 1142 686 1113 780
26 Detox	5000 3176 2489 1552 7344 2842 664
27 Energy Cleansing	7240 7301 run with pads face up in room
28 Epstein Barr Support	105 172 253 274 660 667 669 738 825
29 Fungus/Mold	4442 2411 1833 1823 1333 1155 1130
30 Lymph	432 880 727 3176 478 337 833

This list is NOT claiming to treat, cure, or prevent any disease. It is only meant to provide a frequency environment that supports the body into coming back into natural balance. Please seek proper medical care when needed.



Harmonic Light in Our Energy Field Quantum Effect



Distance Session Client Feedback

“Her BE LIGHT work is phenomenal and she works on them with me virtually as I am in New York. This past week I had a full week of peace and love; no anxiety.”



“Everything is
energy
& that’s all there is to it.

Match the frequency of the
reality you want and you cannot
help but get that reality. It can
be no other way. This is not
philosophy. This is physics. ”

— *Darryl Anka*

BAGONGPINAY
newfilipina.com

Experience More Harmonic Light

- **Sessions**
- **Rental**
- **Purchase**



Sessions




- **Long Distance Quantum Sessions**
- **Personalized Custom Sessions**
- **Typical Cost is \$40-120 per Session.**

Our Trial Program

Experience the benefit of daily sessions of Harmonic Light. We provide education and support. It's as easy as 1-2-3.

- * **1. \$200 Go3/Wands. \$250 Flex 3, \$300 Grow6 or \$350 Pro6**
Use daily for two weeks.
Enjoy the lights several times each day
Share them with your family, friends, or clients.
Typical Cost for Household of 2 is < \$6 per session
- * **2. Free training on how to use the system.**
Exclusive daily email tutorial
Personalized consultations to personalize your program.
- * **3. Apply the trial cost to purchase.**



Choosing A Light System

- FDA OTC Clearance for Pain Relief, Increased Circulation, & Injury Recovery
- Medical Grade Diodes. The more diodes the stronger the transdermal effect.
- Light Delivered is Pulsed at Highly Accurate Harmonic Frequencies with Custom Program Capabilities
- Designed For Multi-Dimensional Health of Mind, Body, and Spirit
- Customer Service with Free, Ongoing Educational Support
- 3 Year Warranty Minimum
- Extended Warranties Even Better – 2 Year Extension Available for 5 years of total coverage.
- Cost Effective Trial Programs
- Pennies per session for daily use.



LumiVibe Go2/Go3

\$1100 Go2 with one pad

\$1500 Go3 with two pads



Cordless Portable Controller

- 12 Red Diodes for Soft Tissue
- 12 Infrared Diodes for Deep Penetration
- ABC 1-7 Pulsation Settings
- 3 Sweeps, 24 Frequencies
- Rechargeable Lithium Ion Battery
- Perfect for travel and athletes




Pad Accessories:

Portable Red/Blue Eye Mask

- 52 Red and 52 Blue Diodes
- Relaxation & mood support
- Recessed diodes for added comfort

Portable Blues Buster Pad:

- 50 Deep Penetrating Infrared Diodes
- 40 Cleansing & Calming Blue Diodes.



LumiPure Chakra Wand Set



- One Wand with Interchangeable Magnetic Pure Color Light Heads

- Red 660nm Root Chakra
- Orange: 620nm Sacral Chakra
- Yellow: 580nm Solar Plexus
- Green: 525nm Heart Chakra
- Turquoise: 490nm High Heart Chakra
- Sky Blue: 460nm Throat Chakra
- Indigo: 432nm Brow or Third Eye
- Violet: 405nm Crown Chakra

- Pure Wavelength Color Delivery
- Gently Pulses at Harmonic 432hz
- Cord-Free
- Rechargeable Battery
- Perfect for all Energy Healing Modalities
 - Color Therapy
 - Foot Zoning
 - Reiki
 - Chakra Balancing & Clearing



\$999

LumiVibe Flex 3:

3 Port Controller

- Easy ABC 1-7 Settings
- 3 Sweeps
- 25 Frequencies

2 Pads Standard

- Relaxation Eye Mask
- Red/Infrared Body Pad
- Flexible Options
 - Add an Additional Pad or Accessory
 - Upgrade a Pad



LumiVibe Grow 6

6 Port Controller

- Digital Display
- Pause/Play & Timer
- ABC 1-7, 25 Frequencies
- Can Upgrade to Pro6

3 Pads Standard

- Red/IR Body Pad
- Red/IR Local Pad
- Red/Blue Eye Mask
- Add or Upgrade Accessories
 - Blue Local or Body Pad
 - Facial Rejuvenator, Chakra Wands
 - Super T Neuro Pad, Pain Buster



Base Price \$3400

The LumiCeutical Pro6 Controller

Discover the Power of Harmonic Resonant Frequencies Delivered by Light Pulses
 Light at 9999 Frequencies
 88 Customizable Programs
 Amazing Ease of Use, Versatility & Support



Harmonic Light Custom Setting Possibilities

Sample set compiled from the [Electroherbalism Frequency List](#) and wellness professionals

Whole Body Support

- 11 General Power Sweep
- 12 Grounding Balance
- 13 Abundance
- 14 Transformation
- 15 Chakra Balancing
- 16 Instant Vacation
- 17 Master Healing Support
 - ___ 17 Chakra Wand Sweep
 - ___ 17 DNA Repair Support
 - ___ 17 Energize
 - ___ 17 Schumann Wave
 - ___ 17 Ultimate Health
 - ___ Also See #94, 98, 99

Mood and Sleep Support

- 18 Peace & Calming
- 19 Insomnia Support
- 20 Depression /Anxiety Support
 - ___ 20 Emotion Sweep
 - ___ 20 Joy & Happiness
 - ___ 20 Sleep Apnea
 - ___ Also See Sweep #97

Immune & Cleanse Support

- 21 A+ Immune Support
- 22 Allergies
- 23 Autoimmune Support
- 24 Colds/Flu Support
- 25 Chemical Detox
- 26 Detox
- 27 Energy Cleansing - run w/ pads face up
- 28 Epstein Barr Support
- 29 Fungus/Mold
- 30 Lymph
- 31 Heavy Metal
- 32 Herpes/Shingles Support
- 33 Immune Support General
- 34 Immune Support Breast
- 35 Immune Support Lymph
- 36 Immune Support Prostate
- 37 Immune Support Blood
- 38 Inflammation/Infection
- 39 Lyme Support
- 40 Ear Infection/Otitis support
- 41 Papilloma virus support
- 42 Pneumonia Support
- 43 Sinus Support
- 44 Strep Support
- 45 Staph Support
 - ___ 45 Bell's Palsy Support
 - ___ 45 Chronic Fatigue Support
 - ___ 45 Hepatitis C Support
 - ___ 45 Immune Stimulation
 - ___ 45 Macrophage - immune
 - ___ 45 Bronchitis Support
 - ___ 45 Chemical Sensitivity
 - ___ 45 EMF Clean
 - ___ 45 Sore Throat
 - ___ 45 Viral Infection Support
 - ___ See Sweep #94 and #95

Digestive Tract Support

- 46 Acidosis/Alkalize
- 47 Candida
- 48 Colitis/ IBS Support
- 49 Digestion
- 50 Gall Bladder Support
- 51 Parasite Support
- 52 Polyp Support
- 53 Ulcer Support
- 54 Urinary Tract
 - ___ 54 Flukes
 - ___ 54 Hernia Support Masters
 - ___ 54 Hiatal Hernia
 - ___ 54 Hookworm
 - ___ 54 Larynx Support Masters
 - ___ 54 Tapeworm
 - ___ 54 Constipation
 - ___ See Sweeps #94 and 95

Heart & Lung Support

- 55 Circulation
- 56 Lungs
- 57 Artery Support
 - ___ 57 Raynauds Support
 - ___ See Sweep #94

Gland & Hormone Support

- 58 Adrenal /Stress Support
- 59 Breast Fibroid Support
- 60 Endocrine Support
- 61 Endometriosis/Uterine
- 62 Hypothyroid Support
- 63 PM5/Menopause Support
- 64 Prostate Support
 - ___ 64 Pineal Opening
 - ___ 64 Adrenal Balance
 - ___ 64 Hot Flashes
 - ___ 64 Hypothalamus
 - ___ See Sweep #96

Diabetic Support

- 65 Blood Sugar
- 66 Pancreas
- 67 Neuropathy Support

Skin Support

- 68 Acne
- 69 Eczema Support
- 70 Gums Support
- 71 Hives
- 72 Warts
- 73 Rosacea Support
- 74 Psoriasis Support
- 75 Skin Support
- 76 Skin Immune Support
- 77 Wounds/Injuries
 - ___ 77 Healthy Skin
 - ___ 77 Insect Bites
 - ___ 77 Scarring
 - ___ 77 Vitiligo Support

Head & Brain Support

- 78 ADD/ Autism Support
- 79 Headache Support
- 80 MS Support - Masters
- 81 Tinnitus Support
- 82 Parkinson Support
- 83 Cognition Support
 - ___ 83 MS Support 2
 - ___ 83 Vertigo
 - ___ 83 ADD Support

Eye Health Support

- 83 Vision - Sharpen
- 83 Macular Support
- 83 Eye Lens Support
- 83 Eye Nerve Support
- 83 Eye Health Support

Joint, Bone & Muscle

- 84 Arthritis Support
- 85 Pain Support
- 86 Lungs
- 87 Knee/Hip Support
- 88 Bone Trauma
- 89 Sciatica Support
 - ___ 89 Gout Support
 - ___ 89 Fracture Support -Masters
 - ___ 89 Disc Support
 - ___ 89 Sprain
 - ___ 89 Elbow Pain
 - ___ 89 Carpal Tunnel
 - ___ See Sweep #94 and 95

Body Sweating

- 90 Weight Release 1 Support
- 90 Belly Fat Release
- 90 Cellulite Release
- 90 Hypothalamus
- 90 Weight Release


Liver & Kidney

- 91 Liver
- 91 Cirrhosis/Liver Support
- 91 Kidney Stones Support

Ott Settings

- 92 Ott Trace Minerals
- 93 Ott Macro Minerals
- 94 Ott Cleanse, Energy, Rejuvenate 95
- 95 Ott Immune, Inflammation Support
- 96 Ott Hormone Balance
- 97 Ott Relaxation/Clear Thinking
- 98 Ott B Alternative
- 99 Ott New C Alternative
- 99 Anahata Solfeggio
- 99 Ott C Original Alternative
- 99 Ott Harmonic Balance
- 99 Ott Ozon
- 99 Ott Solfeggio

The ___ before a number means that it is available in your software USB but must be manually exchanged for one of the existing modes 12-99 on your controller.



*Offered as support for the body's natural healing response. Not intended to treat, prevent, diagnose, or cure disease process. If you have a medical condition, seek advice from your health care professional.

LumiVibe Pro6:

Our Favorite System

- Pro6 Controller
- Digital display, timer
- ABC 1-7
- 88 Custom Sweeps
- 9999 Frequencies
- Unlimited Possibilities



4 Pads Standard


- Relaxation Eye Mask
- Red/Infrared Body Pad
- Blue/Infrared Body Pad
- Red/Infrared Pain Buster



or

Choose a Free Gift \$700 Value:

- Portable Red/IR Handheld
- Facial Rejuvenator
- 60% off Chakra Wand Set (+\$300)
- 50% off Super T Pad (+\$750)



Base Price \$6210

LumiVibe Pro10:

Our Best Value Save \$1500

- 6 Port Pro Controller
- ABC 1-7
- 88 Custom Sweeps
- 9999 Frequencies
- Unlimited Possibilities

8 Pads Standard

- Eye Mask
- Red/Infrared Body Pad
- Blue/Infrared Body Pad
- Red/Infrared Pain Buster
- Red/Infrared Local Pad
- Blue/Infrared Buster Pad
- Renew Facial Rejuvenator
- Go2 Handheld System – Handheld & Portable Eye Mask or Blue Buster
- Mix and match to upgrade to Super T or Include Chakra Wands




Can
Be
Substituted



Pad Accessories

Eye Mask \$585

- 54 Red
- 54 Blue
- Calming for the Nervous System
- Comes in 2 styles – Portable & Standard



Red Local Pad \$750

- 60 Red Diodes
- 72 Near Infrared Diodes




Blue Local Pad \$750

- 60 Blue Diodes
- 72 Near Infrared Diodes

Red Body Pad \$1050

- 120 Red Diodes
- 144 Near Infrared Diodes



Blue Body Pad \$1050

- 120 Blue Diodes
- 144 Near Infrared Diodes

Pad Accessories

Pain Buster \$550

One of our favorite pads with 40 red and 50 near infrared diodes in our most concentrated array. Wonderful for feet, neck, and small areas that need the power of vibrational light.



Blues Buster \$550

40 blue and 50 near infrared diodes for concentrated calming, collagen and pathogen cleansing support. Comes in a portable version to go with our portable handheld controller for \$530.



Super T Neuro Pad

- *Relieve Pain from Neuropathy, Arthritis, Muscle Tension & Injury*
- *Increase Circulation & Cellular Regeneration*
- *Experience a Relaxed State of Wellbeing*



Exclusive design of over 56 Red and 248 Near Infrared diodes
Comfortably Conforms to

- Head
- Sacrum to Spine
- Shoulder to Tailbone
- Foot
- Perfect for the Whole Body



Can be upgraded into our Flex, Grow and Pro Systems
\$1500



shinewithlight.com



RENEW Facial Rejuvenator

Daily Anti-Aging Light Therapy

Red, Near infrared and Orange wavelengths.

Research has found that applying these wavelengths to the skin stimulates fibroblasts to produce collagen.

FDA over the counter (OTC) cleared for the treatment of fine lines and wrinkles.

Combine with additional pads to support increased circulation, cleansing and cell regeneration throughout the body including the liver, the skin's companion organ.

Can be used with the Flex 3, Grow6 or Pro6 Controller.

\$700

shinewithlight.com




Two Weeks That Can Change Your Life

Try before you buy and experience the benefit of daily sessions of Light Energy. We provide education and support. It's as easy as 1-2-3.

- * **1. \$250 Flex 3, \$300 Grow6 or \$350 Pro6**
Use daily for two weeks.
Enjoy the lights several times each day
Share them with your family, friends, or clients.
Typical Cost for Household of 2 is < \$6 per session

- * **2. Receive free training on how to use the system.**
Exclusive daily email tutorial
Personalized consultations to personalize your program.

- * **3. You can apply your rental cost to purchase.**
Applied to purchase made within 10 days of the end of your rental.





Exclusive Dealers of LumiCeuticals in North America. Join Our Ambassador Team

Personalized Customer Service

Free Education & Leadership Support

Opportunity – Spread Light through Direct Sales Commissions

www.shinewithlight.com

Shirley: shirleyjoffs@gmail.com

Marta: marta@shinewithlight.com

shinewithlight.com

